AMICA-GEN

Generatori RF & MW programmabili
Programmable RF & MW generators

CODICE / CODE DESCRIZIONE / DESCRIPTION
AGN-3.0 Generatore MW a stato solido, programmabile, 2450 MHz/140W, con pompa peristaltica integrata
Solid state, programmable MW generator, 2450 MHz/140W, with built-in peristaltic pump

AGN-H-1.0 Generatore programmabile RF (450kHz/200W@500hm) e MW (2450 MHz/140W), con pompa peristaltica integrata
Programmable RF generator (450kHz/200W@500hm) and solid state MW generator (2450 MHz/140W), with built-in peristaltic pump

AMICA-PROBE

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* For the high power version with short coaxial cable (1.5m) add to code the following suffix: -S1.5

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♦ For the detachable cables version add the following suffix: -DC

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* For the Osteoma Osteoid (00) treatment version add to code the following suffix: -OS

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THE NEW GENERATION OF INTERSTITIAL
PROBES FOR MICROWAVE AND
RADIOFREQUENCY THERMOABLATION

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APPARATUS FOR MICROWAVE AND RADIOFREQUENCY THERMOABLATION

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AMICA GEN is a programmable microwave (MW, 2450 MHz) and radiofrequency (RF, 450 kHz) generator, the only available generator capable to provide the two most widespread thermoablation technologies in the same hardware.

Generator features:

- Energy output: up to 140 W continuous wave at 2450 MHz; up to 200 W at 450 kHz on a 50Ohm load.
- Electrical safety: DC-powered MW and RF amplifiers, low dispersed currents.
- Constant energy delivery monitoring: continuous measurement of forward and reflected MW power. Continuous monitoring of RF power, current and impedance. Calculation of the net MW or RF energy delivered to target tissues. MW and RF amplifiers self-diagnostics. Energy delivery automatically stopped in case of alarm or error.
- Essential and intuitive user interface: just an LCD touchscreen and a knob for browsing through menus and setting the operative parameters.
- "Open" digital architecture: any customization simply requires a software upgrade.
- Interfacability: AMICA-GEN may be interfaced to a PC or other external devices through a dedicated serial (RS232) port.
- Automatic peristaltic pump, integrated in the generator.
- Automatic track ablation: AMICA GEN is the only generator featuring a built-in automatic procedure for track ablation.
- Outstanding portability: AMICA GEN is the lightest (12kg) and most compact (45x38x13cm) thermoablation apparatus on the market.

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MICROWAVE ANTENNAS FOR THERMOABLATION

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AMICA PROBE is a disposable applicator for microwave thermoablation of soft tissues.

Probe features:

- Sharp trocar point for easy penetration into tissues.
- Very good echographic visibility.
- Internal water cooling. Pre-assembled input and output hydraulic lines.
- Embedded memory chip for probe identification. Built-in thermocouple.
- Patented device for reflected microwaves entrapment (mini-choke*): no back heating effects! Enhanced control over the coagulation pattern and size, enhanced treatment safety and repeatability.
- No grounding pads.
- Multiple applicators may be operated simultaneously.
- Interstitial thermocouples available for temperature monitoring within or around the lesion.
- Flexible version available for endocavitary use (AMICA-PROBE FLEX, 13G x 1800mm)

Available versions:

gauge: 11G, 14G, 16G
shaft length: 150, 200, 270 mm
built-in or detachable cables
high power version (100W)

AMICA-PROBE coagulative performance on ex vivo bovine liver at room temperature

Ablation size (Length x Diameter)	Ablation time			
	5 min	10 min	15 min	
Delivered power* *net at the antenna tip	20W	27 x 20 mm	31 x 27 mm	38 x 33 mm
	40W	36 x 27 mm	49 x 36 mm	50 x 42 mm
	60W	48 x 34 mm	54 x 40 mm	61 x 48 mm
	80W	52 x 37 mm	66 x 46 mm	73 x 55 mm

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*Worldwide patent by Italian National Research Council (CNR), exclusively licensed to HS.

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RADIOFREQUENCY ELECTRODES FOR THERMOABLATION

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RF AMICA PROBE is a disposable, interstitial RF electrode for soft tissues thermoablation.

Probe features:

- Straight needle (no hooks).
- Monopolar electrode. Grounding pads included in the disposable kit.
- Sharp trocar point for easy penetration into tissues.
- Very good echographic visibility.
- Internal water cooling. Pre-assembled input and output hydraulic lines.
- Embedded memory chip for probe identification.
- Internal thermocouple to monitor the temperature of the tip of the applicator.
- Interstitial thermocouples available for temperature monitoring within or around the lesion.

Probe models already available:

gauge: 17G
shaft length: 100, 150, 200, 270 mm
exposed tip length: 5, 10, 15, 20, 25, 30, 35 mm

RF AMICA-PROBE coagulative performance on ex vivo bovine liver at room temperature

Ablation size (Length x Diameter) Automatic energy delivery mode, 200W max, 10 minutes	Exposed tip length		
	10 mm	20 mm	30 mm
	13 x 11 mm	30 x 22 mm	41 x 33 mm