



## New TL-Rapid Struts Provide Quick Solutions

- Quicker and easier frame assembly
- Acute adjustment with micro-compression or distraction
- Lockable universal hinge on both ends of strut




*TrueLok is Simple  
TrueLok is Stable  
TrueLok is Versatile*



Short Range	Medium Range	Long Range
70 - 104mm	98 - 188mm	163 - 300mm

Struts are available in three sizes:	
Length	Part #
SHORT	50-10170
MEDIUM	50-10180
LONG	50-10190



Manufactured by:  
 ORTHOFIX Srl  
 Via Delle Nazioni 9  
 37012 Bussolengo (Verona)  
 Italy

Telephone +39 045 6719000  
 Fax +39 045 6719380



Your Distributor is: