


# Endovascular Kit

The EndoSet is used for carrying out endovascular laser therapy.


The EndoSet and its individual components (EndoFiber and Insertion Set) are delivered in sterile condition and are meant for single use only.

## EndoFiber (NIRQ600)



1. Fiber, 600 $\mu$ m diameter
2. Clamping adapter
3. Catheter 5fr
4. Ring mark  
(for control by ultrasound)
5. Scale

## Insertion Set



6. Guide wire, 0.35"
7. Insertion lock 5Fr
8. Dilator
9. Angiography needle

Jena Surgical GmbH  
Brüsseler Str. 10  
07747 Jena - Germany  
Phone: +49 (0) 3641 7700 122  
Fax: +49 (0) 3641 7700 102  
Email: [info@jenasurgical.com](mailto:info@jenasurgical.com)  
[www.jenasurgical.com](http://www.jenasurgical.com)

**JENA SURGICAL**  
LASER AT YOUR SIDE