

Components

The Galaxy Fixation components are MR conditional when marked with "▲".

Rods Diameter 12mm ▲

Code	Description
932100	Rod 100mm long
932150	Rod 150mm long
932200	Rod 200mm long
932250	Rod 250mm long
932300	Rod 300mm long
932350	Rod 350mm long
932400	Rod 400mm long

Rods Diameter 9mm ▲

Code	Description
939100	Rod 100mm long
939150	Rod 150mm long
939200	Rod 200mm long
939250	Rod 250mm long
939300	Rod 300mm long

Rods Diameter 6mm ▲



Code	Description
936060	Rod 60mm long
936080	Rod 80mm long
936100	Rod 100mm long
936120	Rod 120mm long
936140	Rod 140mm long
936150	Rod 150mm long
936160	Rod 160mm long
936180	Rod 180mm long
936200	Rod 200mm long

L Rod Diameter 6mm

Code	Description
936010	6mm L Rod

Note: All components of Galaxy Fixation System frames must be identified as MR Conditional prior to being placed in or near an MR Environment. The Threaded Wires (93100), the Wire Locking Clamps (93620), and the L Rod (936010) and Semi-Circular Rods (939010, 939020, 939030) are not MR Conditional. Any construct/frame that is using Threaded Wires, the Wire Locking Clamps, the L-Rod and Semi-Circular Rods must therefore be considered as MR Unsafe.

Clamps



Modules



THE SYSTEM IS MRI CONDITIONAL at 1.5 and 3 Tesla
More details are available on relevant Instructions For Use - PQGAL

Get it all sterile!

Besides the sterile pre-packaged kits, Galaxy Fixation System offers all **clamps and rods single-packed and sterile**, except for the 6mm L-rod cod. 936010. They can be ordered using the item code number preceded by 99- (e.g. 99-939300).

99-93501 Pelvis Sterile Kit

Consisting of:

Code	Description
8x93010	Large Clamp
1x932350	Rod Ø 12mm L 350mm
1x932300	Rod Ø 12mm L 300mm
2x932200	Rod Ø 12mm L 200mm
4x912640	Self-drilling XCaliber Screws, L 260mm, thread length 40mm
1x11138	Drill Guide Ø 4.8mm
1x11137	Screw Guide 80mm
1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect
1x91150	Bone Screw T Wrench

99-93502 Lower Limb Diaphyseal Sterile Kit

Consisting of:

Code	Description
6x93010	Large Clamp
1x932300	Rod Ø 12mm L 300mm
2x932150	Rod Ø 12mm L 150mm
2x912650	Self-drilling XCaliber Screws, L 260mm, thread length 50mm
2x911550	Self-drilling XCaliber Screws, L 150mm, thread length 50mm
1x11138	Drill Guide Ø 4.8mm
1x11137	Screw Guide 80mm
1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect
1x91150	Bone Screw T Wrench

99-93503 Ankle Sterile Kit


Consisting of:

Code	Description
6x93010	Large Clamp
1x932300	Rod Ø 12mm L 300mm
1x932200	Rod Ø 12mm L 200mm
2x932150	Rod Ø 12mm L 150mm
2x911540	Self-drilling XCaliber Screws, L 150mm, thread length 40mm
1x11138	Drill Guide Ø 4.8mm
1x11137	Screw Guide 80mm
1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect
1x91150	Bone Screw T Wrench
1x93080	Transfix Pin 80mm

99-93601 Galaxy Wrist Sterile Kit

Consisting of:

Code	Description
2x93330	Small Multiscrew Clamp-Short
1x93350	Wrist Module
2x19995	Screw Guide
2x13530	Drill guide Ø2.7mm
1x91017	Allen Wrench
4x947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø4mm Thread 3mm L 100/20 QC
2x1-1355001	Drill Bit Ø2,7mm
1x93160	Screw T Wrench QC
3x13715	Kwire Ø1,5x150mm



In a world of trauma
a **GALAXY** of solutions

Orthofix ABS – A touch of Orthofix Knowledge
We are now covering all the mobile devices:
iOS version (App Store) for iPad
Android version (Google Play Store) for Android tablet
Website: <http://abs.orthofix.it/home.asp> for smartphone and PC



Try the app and then insert your review!
This will help us to always innovate.

Try the app and then insert your review!
This will help us to always innovate.

ORTHOFIX Srl
Via Delle Nazioni 9
37012 Bussolengo (Verona)
Italy

Telephone +39 045 6719000
Fax +39 045 6719380

www.orthofix.com
GF-1202-PL-E0 E 02/14




Your Distributor is:




Simple components,


Infinite solutions



WRIST




N.	Code	Description
5	93310	Small Clamp
1	936140	Rod d 6mm L 140mm
1	936080	Rod d 6mm L 80mm
4	36104	Wrist Self-Drilling Cortical Screws Shaft Ø 4mm Thread Ø 3.3-3.0mm 70/20



N.	Code	Description
6	93310	Small Clamp
1	936100	Rod d 6mm L 100mm
1	936080	Rod d 6mm L 80mm
1	936060	Rod d 6mm L 60mm
4	36104	Wrist Self-Drilling Cortical Screws Shaft Ø 4mm Thread Ø 3.3-3.0mm 70/20

WRIST




N.	Code	Description
1	93350	Wrist Module
2	93330	Small Multiscrew Clamp-Short
4	35100	Bone Screws (Shaft Ø 4mm - Thread Ø 3.3-3.0mm)
or		
4	M321	Self-drilling Bone Screws Shaft Ø 3mm - Thread Ø 3.0-2.5mm



N.	Code	Description
1	936010	6mm L Rod
1	93330	Small Multiscrew Clamp-Short
1	93320	Small Multiscrew Clamp-Long
4	35100	Bone Screws (Shaft Ø 4mm - Thread Ø 3.3-3.0mm)
or		
4	M321	Self-drilling Bone Screws Shaft Ø 3mm - Thread Ø 3.0-2.5mm

FEMUR



N.	Code	Description
6	93010	Large Clamp
1	932100	Rod d 12mm L 100mm
2	932200	Rod d 12mm L 200mm
4	912640	XCaliber Screws 260/40 6.0-5.6mm




N.	Code	Description
6	93010	Large Clamp
2	932250	Rod d 12mm L 250mm
1	932350	Rod d 12mm L 350mm
4	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50

FOOT & ANKLE



N.	Code	Description
8	93010	Large Clamp
1	932200	Rod d 12mm L 200mm
1	932100	Rod d 12mm L 100mm
1	932250	Rod d 12mm L 250mm
1	932400	Rod d 12mm L 400mm
1	911540	XCaliber Screw 150/40 6-5.6mm
2	911530	XCaliber Screw 150/30 6-5.6mm




N.	Code	Description
4	93010	Large Clamp
2	99-93030	Large-Medium Transition Clamp (Sterile)
1	932350	Rod d 12mm L 350mm
1	932250	Rod d 12mm L 250mm
1	932200	Rod d 12mm L 200mm
1	92080	Transfix Pin 80mm, Ø 4mm
2	911540	XCaliber Screw 150/40 6-5.6mm

WRIST



N.	Code	Description
3	99-93100	Threaded Wire Sterile (pack of 2)
3	93620	Wire Locking Clamp
4	93310	Small Clamp
1	936080	Rod d 6mm L 80mm
1	936140	Rod d 6mm L 140mm

HUMERUS




N.	Code	Description
6	93110	Medium Clamp
2	939100	Rod d 9mm L 100mm
1	939150	Rod d 9mm L 150mm
4	911530	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/30

WRIST




N.	Code	Description
1	93020	Multiscrew Clamp
1	932200	Rod d 12mm L 200mm
1	93410	Elbow Hinge
1	939150	Rod d 9mm L 150mm
2	93110	Medium Clamp
2	911530	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/30
2	10137	Cortical Screws Shaft Ø 6mm Thread Ø 4.5-3.5mm 120/20

ELBOW




N.	Code	Description
5	93110	Medium Clamp
1	939100	Rod d 9mm L 100mm
1	939150	Rod d 9mm L 150mm
4	10137	Cortical Screws Shaft Ø 6mm Thread Ø 4.5-3.5mm 120/20

WRIST



N.	Code	Description
9	93010	Large Clamp
3	932300	Rod d 12mm L 300mm
2	932250	Rod d 12mm L 250mm
4	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50

PELVIS



N.	Code	Description
8	93010	Large Clamp
2	932200	Rod d 12mm L 200mm
1	932300	Rod d 12mm L 300mm
1	932400	Rod d 12mm L 400mm
4	912650	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 260/50

WRIST



N.	Code	Description
6	93010	Large Clamp
1	932100	Rod d 12mm L 100mm
1	932150	Rod d 12mm L 150mm
1	932200	Rod d 12mm L 200mm
4	911540	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/40

TIBIA & KNEE



N.	Code	Description
8	93010	Large Clamp
2	932350	Rod d 12mm L 350mm
1	932250	Rod d 12mm L 250mm
4	911540	XCaliber Screw 150/40 6-5.6mm
2	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50

THE SYSTEM IS MRI CONDITIONAL AT 1.5 AND 3 TESLA