



Proctoscope (19 mm / 25 mm)

The Evexar Self-Illuminating 19 mm and 25 mm Proctoscopes are single use, fully disposable devices that utilise an illuminator and a light rod to direct light to the end of the scope where it is most needed. This enables examinations and minor therapeutic procedures to be performed in all areas such as out-patients' departments, Intensive Care, Accident and Emergency, or GP surgeries, as well as in the operating room, giving complete clinical freedom to use anywhere.

- **Ideal for Specialist Clinicians for use in Hospitals and Private Clinics, as well as for GP's in their surgeries**
- **Self-illuminated; no need for fibre optic cables or cleaning**
- **Integrated obturator design prevents obstruction of the examination channel and light guide (patent pending)**
- **Fully disposable as a unit after use**
- **Reduced risk of cross-infection since the scope is single use and disposable**
- **Convenient and cost effective**

Order No.	Sales Unit
19 mm: 31076	25 pcs
25 mm: 31077	25 pcs

