

Instrument Set for Removal of Broken Screws **59-050.00**



Aluminium Case, small, silver including tray for instruments (2 pcs.)



Hollow Reamers complete to expose broken screws (cutting counter-clockwise)

- 1 x for 3.5 mm dia. Screws
- 1 x for 4.5 mm dia. Screws
- 1 x for 6.5 mm dia. Screws



Extracting Bolts to remove broken screws (with inner counter-clockwise thread)

- 1 x for 3.5 mm dia. Screws
- 1 x for 4.5 mm dia. Screws
- 1 x for 6.5 mm dia. Screws



Conical Extraction Screw to remove with damaged hexagonal sockets (with counter-clockwise threads)

- 1 x for 3.5 mm dia. Screws, hexagonal recess
- 1 x for 4.5 mm dia. Screws, hexagonal recess
- 1 x for 6.5 mm dia. Screws, hexagonal recess

General Instruments

- 1 x Hollow Gouge to expose broken screws
- 1 x Forceps for broken Screw removal 200 mm
- 1 x Tap Handle with quick coupling 160 mm
- 1 x Sharp Hook
- 1 x Aluminium Template manual for use