

Vacuum System with Plaster Saw
95-460.25

including:

- Vacuum Cleaner + Dust Removal Unit
- Tube, Hose clamp and Tools
- Plaster Saw 95-460.22

not including:

- Magnetic Plate (can be ordered separately)
- Hanging kit for Saw (can be ordered separately)



Main connection	220-240 Volt 50/60 Hz or 110-120 Volt 50/60 Hz
Rated current consumption	1300 Watt
Vacuum	244 mbar / 24,4 KPa
Container Volume	15 l
Filter Performance	99,97% > 0,5 µm
Noise Level	59 dB (A)
Weight	8,0 kg
Cable Length	6 m
Air quantity	61 l/s
Protection Class	II
Safety IP Code	IP2X
Warranty	18 months

Hose clamp

Dust removal unit

Magnetic Plate

Hanging kit for Saw

