


**STERILE CONTAINER SYSTEMS**




**Window Model**

**» 3/4 Quarter Size**

### 3/4 WINDOW MODEL




Specifications		Lid (Perforated)
		Box (Perforated)
h = 100 mm		K6-240D
		K6-241D
		K6-242D
		K6-243D
		K6-244D
		K6-245D
h = 135 mm		K6-246D
		K6-247D
		K6-248D
		K6-249D
		K6-250D
		K6-251D
h = 150 mm		K6-252D
		K6-253D
		K6-254D
		K6-255D
		K6-256D
		K6-257D



Ps: All container measurements cover the outer size of the lids.

**BOX** (perforated)



➔ Should be used with K2-1522 wire base.

code	outside dimension	inside dimension	weight
K6-300D	485 x 272 h=115 mm	422 x 258 h=66 mm	2.00 kg
K6-301D	485 x 272 h=131 mm	422 x 258 h=84 mm	2.16 kg
K6-302D	485 x 272 h=153 mm	422 x 258 h=105 mm	2.26 kg

**PERFORATED LIDS**



Dimension= 477 x 285 mm

code	color	weight
K6-340D		1.30 kg
K6-341D		
K6-342D		
K6-343D		
K6-344D		
K6-345D		

STANDARD MODEL 1

BIO-BARRIER MODEL 33

Polymers PPSU Lids 65

WINDOW MODEL 85

B<sub>3</sub> Sterile MODEL 103

for STERRAD 135

Dental, Mini, Special 145


BASKETS 163

Autolock Mechanism Emulating Indicator 183

ACCESSORIES 191

BASKET MADE FROM STAINLESS STEEL WIRE MESH


WIRE BASKET



code	outside dimension
K2-1500	405 x 250 x h=30 mm
K2-1501	405 x 250 x h=50 mm
K2-1502	405 x 250 x h=70 mm
K2-1503	405 x 250 x h=100 mm



WIRE BASKET WITH LID




code	outside dimension
K3-1030	405 x 250 h=35 mm
K3-1031	405 x 250 h=55 mm
K3-1032	405 x 250 h=75 mm
K3-1033	405 x 250 h=105 mm
K3-1034	ONLY LID, for baskets




WIRE BASKET WITH PERFORATED METAL PLATE

WIRE BASKET



code	outside dimension
K2-1501 P	405 x 250 x h=50 mm
K2-1502 P	405 x 250 x h=70 mm
K2-1503 P	405 x 250 x h=100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1031 P	405 x 250 x h=55 mm
K3-1032 P	405 x 250 x h=75 mm
K3-1033 P	405 x 250 x h=105 mm
K3-1034 P	ONLY LID, for baskets







## BASKET MADE FROM STAINLESS STEEL SHEET

### BASKET (with stand)



	code	outside dimension
 Box and Lid	K3-728	405 x 253 h=43 mm
	K3-729	405 x 253 h=63 mm
	K3-730	405 x 253 h=83 mm
	K3-731	405 x 253 h=113 mm
 Box	K3-748	405 x 253 h=40 mm
	K3-749	405 x 253 h=60 mm
	K3-750	405 x 253 h=80 mm
	K3-751	405 x 253 h=110 mm
 Lid	K3-762	ONLY LID, for baskets



### STAND



code	height
K3-767	20 mm



### STAND



code	height
K3-768	10 mm



## IDENTIFICATION (aluminium)

### LABEL HOLDER



code	color
K2-560	Red
K2-561	Blue
K2-562	Green
K2-563	Yellow
K2-564	Grey
K2-565	Black

### LABELS



code	color
K2-570	Red
K2-571	Blue
K2-572	Green
K2-573	Yellow
K2-574	Grey
K2-575	Black

Coloured identification labels  
Engraved to your specification



⦿ Grey colored recommended for use in Sterrad.  
Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

STANDARD MODEL	1
BIO-BARRIER MODEL	33
Polymers PPSU Lids	65
WINDOW MODEL	85
B <sub>3</sub> Barrier MODEL	103
for STERRAD	135
Dental, Mini, Special	145
BASKETS	163
Autolock Mechanism Emulating Indicator	183
ACCESSORIES	191

**DIVIDERS**



K2-1540  
405 x 250 x 80 mm  
K2-1541  
405 x 250 x 110 mm  
K2-1542  
405 x 250 x 130 mm

**WIRE BASE**



K2-1522  
405 x 255 x 30 mm

**DIVIDERS**



K2-1530  
405 x 250 x 80 mm  
K2-1531  
405 x 250 x 110 mm  
K2-1532  
405 x 250 x 130 mm

(-)

**DRAPE RETAINER**



K2-1521  
405 x 255 x 30 mm

**DRAPE RETAINER**



K2-1520  
405 x 255 mm

**ALUMINIUM BASKET (with lid)**



K3-1119  
245 x 130 x 55 mm

**ALUMINIUM BASKET (with lid)**




K3-1120  
245 x 130 x 70 mm

**ALUMINIUM BASKET (with lid)**




K3-1121  
245 x 130 x 105 mm

**SILICONE MAT**



K2-524  
380 x 230 mm  
K2-525  
440 x 230 mm







**SILICONE MESH MAT**



K2-531  
480 x 250 mm

**SILICONE HANDLE**



K2-710   
K2-711   
K2-712   
K2-713   
K2-714   
K2-715 

**PAPER FILTER**



K2-740  
Single use 100 pcs  
with indicator. Aygun mark.  
• Diameter: 185 mm

**TEFLON FILTER**



K2-742  
Number of sterilization  
cycle > 2000  
• Diameter: 194 mm







**FILTER DISPENSER**




K2-730

**IDENTIFICATION LABELS (aluminium)**




K2-700   
K2-701   
K2-702   
K2-703   
K2-704   
K2-705 

**LABELS (paper)**



K2-673  
100 pcs. With indicator.  
For small size and  
special containers

**LABELS (paper)**



K2-672  
100 pcs. With indicator.  
For full and  
Half size containers


(-)

**AUTO LOCK MECHANISM**



K1-900  
(2 Pcs)

**EMULATING INDICATOR CARD (class 6)**



K2-674  
Paper 250 Pcs.



**SECURITY SEAL (blue)**

K2-663  
100 pcs. with indicator.  
made of plastic



**SECURITY SEAL (blue)**

K2-664  
100 pcs. made of plastic



**SECURITY SEAL (blue)**

K2-661  
100 pcs. with number.  
made of plastic



**SECURITY SEAL (red)**

K2-662  
100 pcs. with indicator.  
made of plastic



**SECURITY SEAL (red)**

K2-660  
100 pcs. made of plastic



**SECURITY SEAL (red)**

K2-665  
100 pcs. with number.  
made of plastic



**DISPOSAL STAMP**

K2-717  
40 x 35 x 9 mm  
Requested number or  
letter must be informed  
during the order.



**LABEL HOLDER**

K2-716  
50 x 35 x 13 mm  
Holder for identification  
labels

(-)