

CODE	DIAMETER	LENGTH	BOX
HS NOTA14100	14G	100mm	10
HS NOTA14150	14G	150mm	10
HS NOTA14200	14G	200mm	10
HS NOTA16100	16G	100mm	10
HS NOTA16150	16G	150mm	10
HS NOTA16200	16G	200mm	10
HS NOTA16250	16G	250mm	10
HS NOTA16300	16G	300mm	10
HS NOTA18100	18G	100mm	10

CODE	DIAMETER	LENGTH	BOX
HS NOTA18150	18G	150mm	10
HS NOTA18200	18G	200mm	10
HS NOTA18250	18G	250mm	10
HS NOTA18300	18G	300mm	10
HS NOTA20100	20G	100mm	10
HS NOTA20150	20G	150mm	10
HS NOTA20200	20G	200mm	10
HS NOTA20250	20G	250mm	10
HS NOTA20300	20G	300mm	10



nota

DISPOSITIVO AUTOMATICO PER
BIOPSIA ISTOLOGICA TRU-CUT

*AUTOMATIC DEVICE FOR TRU-CUT
HISTOLOGICAL BIOPSY*



HS HOSPITAL SERVICE S.P.A

Sede legale/Registered office:
Via Zosimo 13 - 00178 Roma, ITALY

Stabilimento di produzione/Factory:
Via Angela Vacchi 23/25
04011 Aprilia (Latina), ITALY

Tel.: +39 06 9201961
Fax: +39 06 92727871
email: hs@hshospital-service.com
www.hshospital-service.com

NOTA

DISPOSITIVO AUTOMATICO PER BIOPSIA ISTOLOGICA TRU-CUT

AUTOMATIC DEVICE FOR TRU-CUT HISTOLOGICAL BIOPSY

NOTA è un ago automatico a ghigliottina per biopsia dei tessuti molli.

Caratteristiche:

- leggera, maneggevole, facile da usare e da ricaricare.
- caricamento a rotazione.
- indicatore del caricamento.
- campionamento automatico.
- dispositivo di sicurezza attivo che si inserisce automaticamente durante la carica.
- pulsante di sparo posteriore.
- avanzamento 22mm.
- marker ecogeno interno.


NOTA is an automatic soft-tissues tru-cut biopsy needle.

Features:

- light, easy handling, easy to use and to charge.
- loading with rotating system.
- loading indicator.
- automatic sampling.
- security feature that automatically activates during the device charging.
- rear firing position.
- 22mm penetration depth.
- inner echogenic marker.



LEVA DI CARICAMENTO
TRIGGERING LEVER



FERMO DI SICUREZZA
SAFETY BLOCK



AGO CON PUNTA TRACIANTE
CUTTING TIP



TRU-CUT Type

