


# TEARPROOF TISSUE RETRIEVAL BAG



Instrument before use

## Tearproof Retriever Pouch for collecting excised tissue in minimal invasive surgery

The Tissue Retrieval Bag is comprised of a flexible plastic bag with a large, easily-accessible opening, a push/pull rod with thumb ring handle, finger rings, string and suture and an introducer tube. In the fully deployed condition, the opening is maintained in a fully-open position by a metallic rim. A string with a future facilitates closure of the specimen bag after the specimen has been collected.

Order No.	Bag Diameter × Length	Diameter Introducer Tube	Volume (ml)	Sales Unit
QW-3	50 mm × 100 mm	10 mm	140	10
QW-4	50 mm × 150 mm	10 mm	250	10
QW-5	60 mm × 200 mm	10 mm	320	10
QW-6	90 mm × 240 mm	12 mm	1200	10

CE<sub>0197</sub>