


HEMORRHOIDS STAPLER

Procedure for Prolapse and Hemorrhoids (PPH)



THE NEW GENERATION


Innovations → Advantages for the user

- No overcompression of tissue possible → lower tissue trauma
- Easier locking → for better control
- Loud click sound → acoustic signal indicating the correct anastomozation
- Anti-slip covering on anal dilatator → better fixing for a more equal purse-string suture
- Bigger cartridge lumen → for bigger specimen resection
- Transparent anoscope → better control/view when applying the purse-string suture

For S.T.A.R.R. Procedure

The instrument with the additional labelling S.T.A.R.R. contains staples of a length of 4.5 mm for tissues with a thickness of 2.0 mm. Consequently, it is appropriate for complete resections of all tissue layers.

Enclosed Surgical Accessories



Article	Order No.	Number of Staples	Outer Diameter*	Inner Diameter*	Leg Height*	Sales Unit
GH-Y-32	GH-Y-32	30	32	24,2	3,8	3 pcs
GH-Y-34	GH-Y-34	32	34	26,2	3,8	3 pcs
GH-Y-36	GH-Y-36	34	36	28,2	3,8	3 pcs
S.T.A.R.R. Stapler						
S.T.A.R.R.-System GH-34	GH-ST-34	32	34	26,2	4,5	3 pcs
S.T.A.R.R.-System GH-36	GH-ST-36	34	36	28,2	4,5	3 pcs
Surgical Accessories (content see above)						
Hemorrhoid Stapler Set	GH-AS-01	for GH-Y-32				3 pcs
Hemorrhoid Stapler Set	GH-AS-02	for GH-Y-34				3 pcs
Hemorrhoid Stapler Set	GH-AS-03	for GH-Y-36				3 pcs

CE0481

*mm