



# CIRCULAR STAPLER




Tangible turning marks reduce the chance of a mal-operation.




Unique hardened chopping ring provides the precise cutting of the tissue.




Safe and solid closure of anvil and cartridge.




Visible Indicator makes it easier for the operator to keep the triggering range between 1.5 mm – 2.5 mm.



Prevention of a stenosis of the anastomosis by a larger inner lumen.



For different calibers the corresponding staple legs will be triggered into the right figure of »B« for the different tissue thicknesses.



| Article | Order No. | Colour of Protector | Number of Staples | Outer Diameter* | Inner Diameter* | Leg Height* | Sales Unit |
|---------|-----------|---------------------|-------------------|-----------------|-----------------|-------------|------------|
| WH-Y-17 | WH-Y-17   | white               | 16                | 17.0            | 10.0            | 3,8         | 3 pcs      |
| WH-Y-21 | WH-Y-21   | yellow              | 20                | 21.0            | 14.0            | 4,5         | 3 pcs      |
| WH-Y-24 | WH-Y-24   | grey                | 18                | 24.0            | 15.0            | 4,5         | 3 pcs      |
| WH-Y-26 | WH-Y-26   | red                 | 20                | 26.0            | 17.0            | 4,5         | 3 pcs      |
| WH-Y-29 | WH-Y-29   | blue                | 24                | 29.0            | 20.0            | 4,8         | 3 pcs      |
| WH-Y-32 | WH-Y-32   | green               | 26                | 32.0            | 22.0            | 5,0         | 3 pcs      |
| WH-Y-34 | WH-Y-34   | purple              | 28                | 34.0            | 24.0            | 5,0         | 3 pcs      |